

Evaporative Light Scattering Detectors from Polymer Laboratories



Introduction

- Polymer Laboratories has over 15 years of experience in ELSD
- ELSD can outperform traditional detectors when analysing non-chromophoric samples by HPLC
- Traditional HPLC detectors such as UV and RI have limited capabilities:
 - UV and RI are not compatible with a wide range of solvents
 - RI detection is not gradient compatible
 - Different analytes produce different UV responses depending on their extinction co-efficient
- ELSDs can detect anything that is less volatile than the mobile phase
- ELSD is universal and compatible with a wide range of solvents



Introduction

Polymer Laboratories presents 3 models of ELSD, each offering high sensitivity for a wide range of applications:

PL-ELS 2100

- Regular HPLC with low temperature operation for volatile compounds

PL-ELS 1000

- Routine HPLC and GPC using high boiling point solvents

PL-ELS 1000 μ

- For microbore and capillary LC






New PL-ELS 2100 Detector



Principles of Operation


The ELSD principle of operation employs three distinct stages:

- Nebulisation
- Evaporation
- Detection




Nebulisation

- Eluent flow mixed with N₂ or Air
- Concentric nebuliser
- Efficient nebulisation: stable droplet plume, uniform droplet size
- Temperature independently controlled
Narrower cone of the plume
(minimise band broadening)
- Small nebuliser chamber
(reduces band broadening)




Evaporation

- Droplets pass through heated drift tube
 - 30cm straight tube
- Removes mobile phase to leave particulate form of analyte
- Temperature set according to analyte volatility
- Temperature controlled by user
- Important to have laminar flow (reduces band broadening)



Detection

- Particles irradiated with light source
 - LED (ca 400-500 nm)
- Particles scatter light according to their size (mass sensitive)
- Scattered light is detected by photomultiplier at fixed angle from incident light
- Detection independent of optical properties of analyte



Advantages of Evaporative Light Scattering Detection

- Universal - responds to all compounds in the mobile phase
- Not dependent on spectroscopic properties of analyte
- Produces more uniform detection sensitivity for analytes
- Not susceptible to baseline drift during gradient elution, temperature or solvent pump fluctuations
- ELSD compatible with a much wider range of solvents compared to Refractive Index



Advantages of Evaporative Light Scattering Detection

- Removes the need for derivatisation steps (eg amino acids, toxins)
- Fast setup and equilibration
- Sensitivity in the order of 1-50ng (**on column**)
(depends on eluent flow rate)
- No interference from solvent front peaks
(enables fast analysis)
- Removal of mobile phase eluent allows rapid HPLC gradients
- Flow rates up to up to 5ml/min can be achieved with no affect on baseline stability
- Ideal for High Throughput Screening



Fast Gradient, Fast Flow Rate Capability

Column: PLRP-S 100Å 5µm, 50x4.6mm
Eluent A: Water + 0.05% TFA
Eluent B: ACN + 0.05% TFA
Gradient: 5-95% B in 1 min
Detector: PL-ELS 2100
(neb=30°C, evap=30°C, gas=1.6 SLM)


Sample: Indapamide(IND), Dibutyl phthalate (DBP)

Flow rate increased from 2ml/min up to 5ml/min

Note: IND is non-volatile, DBP is relatively volatile




Fast Gradient, Fast Flow Rate Capability




Fast Gradient, Fast Flow Rate Capability

Stable baseline through the gradient



High Throughput Screening

PL-ELS 2100 is transparent to DMSO at ambient temperature !!



Nebuliser 30°C Evaporator 30°C
(peak 2 DMSO)



Evaporative Light Scattering Detection Ideal Complement to LC-MS

- Similar operating principles to LC-MS
 - Volatile buffers
 - Favours lower flow rates (ie 0.2-0.5ml/min)
- Can develop LC methods on ELSD then transfer to LC-MS
- ELSD can provide supporting information when used in tandem with LC-MS



Ideal Complement to LC-MS

Sample Mixture of known 1:1 ratio

LC-MS
results show ratio to be
3:1

UV-Vis
results show ratio to be
10:1


PL-ELS 2100
results show ratio to be
1:1
(Response independent
of optical properties)

Sample Report:


Sample 1 Vial 1:

Date 14-Jul-2003 Time 11:25:22 Description Testm


1: MS ES+ :TIC Smooth (Mn, 2x2)



UV



ELSD



Operation of the PL-ELS 2100


For Volatile & Semi-Volatile Solutes

- A new approach to ELSD applications
- Nebuliser and evaporator temperatures determined by the volatility of the ANALYTE
- Gas flow used to aid evaporation of the mobile phase
 - 100% water @30°C requires high gas flow
 - 100% Hexane @30°C requires no gas flow
- PL-ELS 2100 can operate in 100% water @ 25°C



Non-Volatile Application: Sugars

Column: PL Hi-Plex Ca 9 μ m, 300x7.7mm
Eluent: Water
Flow Rate: 0.6ml/min
Temp: 85°C
Inj Vol: 10 μ l
Detector: PL-ELS 2100 (neb=30°C, evap=90°C, gas=1.6 SLM)
Samples: 1. Fructose, 2. Glucose, 3. Sucrose I, 4. Lactose, 5. Stachyose



Non-Volatile Application: Phospholipid Separation

Column: Lichrospher DIOL 5 μ m,150x2.1mm
Eluent A: IPA/Hexane/Water/Ammonia Hydroxide
57.8/40/2/0.2


Eluent B: IPA/Hexane/Water/Ammonia Hydroxide
51.8/40/8/0.2

Gradient: 0-100% B in 7 mins (hold 8 mins)
100-0% B in 5 mins (hold 10mins)

Flow Rate: 0.3ml/min

Detector: PL-ELS 2100
(neb=30°C, evap=80°C, gas=1.0 SLM)

Samples: 1. Cholesterol
2. Phosphatidylethanolamine
3. Phosphatidylcholine
4. Sphingomyelin
5. Lysophosphatidylcholine




Semi-Volatile Application: Underivatised Amino Acids

Column: Thermo-Hypersil ODS 5 μ m, 250x4.6mm
Eluent A: Water
Eluent B: Acetonitrile
Gradient: 100% A in 5 mins hold, 0-40% B in 20 mins
Flow Rate: 0.6ml/min
Inj Vol: 10 μ l
Detector: PL-ELS 2100
(neb=50°C, evap=50°C, gas=1.6 SLM)



Semi-Volatile Application: Underivatised Amino Acids

ELSD removes the need for derivatisation for applications such as amino acids



Volatile Application: Effect of Changing Temperature

Column: Adsorbosil C18 5 μ m, 150x4.6mm
Eluent A: Water + 0.1% TFA
Eluent B: ACN + 0.1% TFA
Gradient: 60-90% B in 5 mins
Flow Rate: 1.0ml/min
Detector: PL-ELS 2100
(neb/evap=same temperature, gas=1.8 SLM)

Samples:


1. Acetanilide
2. Indapamide
3. Ibuprofen
4. Dibutylphthalate

Note: Peak 2 is non-volatile
Peaks 1, 3 & 4 are relatively volatile



Volatile Application: Effect of Changing Temperature

- At high temperature, peak 2 has better S/N but peaks 1,3 & 4 are not detected
- Running cold, all four peaks are detected




Volatile Application: Phthalate Separation

Column: X-Terra C18 3.5 μ m, 30x2.1mm
Eluent A: Water
Eluent B: ACN
Gradient: 0-100% B in 3 mins (hold 0.7 min)
Flow Rate: 0.5ml/min
Detector: PL-ELS 2100
(neb=25°C, evap=25°C, gas=1.6 SLM)

Samples


1. Diethylphthalate
2. Dipropylphthalate
3. Dibutylphthalate
4. Dipentylphthalate
5. Dioctylphthalate



Volatile Application: Parabens Separation

Column: Adsorbosil C18 5 μ m, 150x4.6mm
Eluent A: Water
Eluent B: ACN
Gradient: 50-75% B in 7 mins
Flow Rate: 1.0ml/min
Detector: PL-ELS 2100
(Evap=30°C, Gas=1.0 SLM)


Samples
1. Methylparaben
2. Ethylparaben
3. Propylparaben



PL-ELS 2100

Detection of Low Molecular Weight Compounds

Column: Adsorbosil C18 5 μ m, 150x4.6mm
Eluent: 60/40 Water/Acetonitrile
Flow Rate: 1.0ml/min
Inj Vol: 10 μ l
Detector: PL-ELS 2100
(neb=30°C, evap=30°C, gas=1.6 SLM)
Sample: Urea 1mg/ml



PL-ELS 2100

Sensitivity: Limits of Detection

Column: Spherisorb C18 5 μ m, 150x4.6mm

Eluent A: 99% Water + 0.1% TFA


Eluent B: 1% ACN + 0.1% TFA

Flow Rate: 0.5ml/min

Detector: PL-ELS 2100

(neb=50°C, evap=70°C, gas=1.2 SLM)


Sample: DL-Glutamic Acid



PL-ELS 2100

Sensitivity: Limits of Detection

Column: ACE C8 3 μ m, 50x4.6mm
Eluent A: 100% Water
Eluent B: 100% ACN
Isocratic: 50/50 A/B
Flow Rate: 1.0ml/min
Detector: PL-ELS 2100
(neb=50°C, evap=50°C, gas=1.4 SLM)
Sample: Caffeine



Additional Benefits of the PL-ELS 2100

- Low temperature operation even for 100% water
- Extremely low dispersion for high resolution separations
- High eluent flow rates, up to 5ml/min
- Improved uniformity of response across a solvent gradient
- Rapid equilibration
- Extremely small footprint, stackable
- Ergonomic design
- Easy to use



Les produits de nos partenaires

Systeme LC ECD ALexys



Système LC ECD ALEXYS



Système intégré dédié à la détection électrochimique:

- Flexibilité : existe en version LC ou μ LC
- Sensibilité:
- Ensemble du système optimisé pour la détection électrochimique
- Compatibilité CFR 21 Par 11

Composants du système :

- Pompe, amortisseur de pulsations
- Dégazeur
- Injecteur automatique analytique ou micro
- **DECADE II**
- Système de traitement de données



Decade II

Véritable plate forme de travail

ANTEC
LEYDEN


- Un four, Cage de faraday, abritant cellule de détection et colonne
- Pouvant contrôler une ou deux cellules de détection (mode série ou parallèle)
- Concept de cellule original (vaste choix d'électrodes de travail et de référence)
- Filtrage efficace du bruit de fond (ADF)



Cellule ampérométrique VT-03

ANTEC
LEYDEN


- Cellule de type Wall Jet
- Stabilisation rapide,
- Faible volume(0,3 μ L à 10nL)
- Sensibilité



A: couche mince, B wall jet, C Poreuse

ADF advanced digital filtration

Filtre conçu spécialement pour la détection électrochimique



Les basses fréquences (pics) passent

Les fréquences élevées (bruit) sont bloquées



Applications types



- Analyse de traces nécessitant: sélectivité et sensibilité
- Applications en Micro dialyse (Off line ou On line)



Micro dialysate (frontal cortex + nucleus accumbens) settings

■ALEXYS system

- flow 1 ml/min; injection 25 ul (ul pick-up mode)
- column 10 x 0.46 cm 3 um ODS-2; spiked with 10 fmol DA + 5HT/25ul (= 4×10^{-10} M)
- ADF 0.01 Hz up to 10 min, thereafter 0.002 Hz

